

TAKTVOLL®

**NEW PRODUCT!** ES-100V PRO LCD  
Touchscreen Electrosurgical System  
**ES-100V Pro**




The ES-100V Pro satisfies the veterinarian's demands with precision, safety, and reliability.



# Vessel Sealing Generator

## ES-100V Pro For Veterinary Use



- Glass full version design, which is extremely tense and beautiful
- Color touch screen operation panel, flexible and easy to operate
- 7 Working modes
- Large blood vessel sealing function
- Electrosurgical Pens and Foot Switch Control
- Recently Mode, Power and Other Parameters
- Quick Adjustment of Power
- Cut and Coag in an Intermittent Manner



### ES-100V Pro Feature

3 Cut Mode		2 Coag Mode		1 Vessel Sealing	
<b>PURE</b>	Provides a standard Cutting effect.	<b>SPRAY</b>	Designed to provide hemostasis in a broader area with less depth of tissue damage when compared to desiccation.	<b>Vessel Sealing</b>	
<b>BLEND 1</b>	Provides a combination of 70% Cutting and 30% Coagulation	<b>FORCED</b>	Fast and effective standard coagulation with moderate thermal penetration. Slight carbonization effects.	Seal and cut vessels up to and including 7mm in diameter. Sealing vessels and tissue bundles without dissection or isolation , easy to operate.	
<b>BLEND 2</b>	Provides a combination of 50% Cutting and 50% Coagulation			<b>1 Bipolar</b>	
				<b>Standard</b>	
				Bipolar coagulation which gently coagulates tissue for minimal tissue damage reducing carbonization and sticking.	